

Product Datasheet

RPA32/RPA2 (Phospho-phospho T21) Conjugated Antibody (orb1612490)

Description	RPA32/RPA2 (Phospho-phospho T21) Conjugated Antibody
Species/Host	Rabbit
Conjugation	Biotin
Storage	Store at 4°C for 6 months.
Note	For research use only
Clonality	Recombinant
Source	Recombinant Rabbit Monoclonal Ab
Uniprot ID	P04629
Dilution Range	Working with enzyme-conjugated streptavidin, 1: 50 - 1: 1,000
Expiration Date	12 months from date of receipt.