

Product Datasheet

P-V-Myb+C-Myb(Phospho-S11) Conjugated Antibody (orb1612534)

Description	P-V-Myb+C-Myb(Phospho-S11) Conjugated Antibody
Species/Host	Rabbit
Conjugation	Biotin
Storage	Store at 4°C for 6 months.
Note	For research use only
Clonality	Recombinant
Source	Recombinant Rabbit Monoclonal Ab
Uniprot ID	Q15661
Dilution Range	Working with enzyme- conjugated streptavidin, 1: 50 - 1: 1,000
Expiration Date	12 months from date of receipt.