

## Product Datasheet

# MYH-pan Polyclonal Conjugated Antibody (orb1618296)

<b>Description</b>	MYH-pan Polyclonal Conjugated Antibody
<b>Species/Host</b>	Rabbit
<b>Conjugation</b>	Biotin
<b>Storage</b>	Store at 4°C for 6 months.
<b>Note</b>	For research use only
<b>Clonality</b>	Recombinant
<b>Source</b>	Rabbit Polyclonal Ab
<b>Uniprot ID</b>	<a href="#">Q9UKX2</a> , <a href="#">P13535</a> , <a href="#">P12883</a> , <a href="#">Q9Y623</a> , <a href="#">P13533</a> , <a href="#">P11055</a> , <a href="#">P12882</a>
<b>Dilution Range</b>	Working with enzyme-conjugated streptavidin, 1: 50 - 1: 1,000
<b>Expiration Date</b>	12 months from date of receipt.