

Product Datasheet

PTTG1/2/3 Polyclonal Conjugated Antibody (orb1618553)

Description	PTTG1/2/3 Polyclonal Conjugated Antibody
Species/Host	Rabbit
Conjugation	Biotin
Storage	Store at 4°C for 6 months.
Note	For research use only
Clonality	Recombinant
Source	Rabbit Polyclonal Ab
Uniprot ID	O95997
Dilution Range	Working with enzyme-conjugated streptavidin, 1: 50 - 1: 1,000
Expiration Date	12 months from date of receipt.