

Product Datasheet

HOXA1/B1/D1 Conjugated Antibody (orb1623884)

Description HOXA1/B1/D1 Conjugated Antibody

Species/Host Rabbit

Conjugation Biotin

Storage Store at 4°C for 6 months.

Note For research use only

Clonality Recombinant

Source Rabbit Polyclonal Ab

Uniprot ID [P49639](#), [Q9GZZ0](#), [P14653](#)

Dilution Range Working with enzyme-conjugated streptavidin, 1: 50 - 1: 1,000

Expiration Date 12 months from date of receipt.