

Product Datasheet

PAK1/2/3 (Ab-423/402/421) Conjugated Antibody (orb1624275)

Description	PAK1/2/3 (Ab-423/402/421) Conjugated Antibody
Species/Host	Rabbit
Conjugation	Biotin
Storage	Store at 4°C for 6 months.
Note	For research use only
Clonality	Recombinant
Source	Rabbit Polyclonal Ab
Uniprot ID	Q13153, O75914, Q13177
Dilution Range	Working with enzyme-conjugated streptavidin, 1: 50 - 1: 1,000
Expiration Date	12 months from date of receipt.