

Product Datasheet

SMC1A(N-term) Conjugated Monoclonal Antibody (orb1629079)

Description	SMC1A(N-term) Conjugated Monoclonal Antibody
Species/Host	Mouse
Conjugation	Biotin
Storage	Store at 4°C for 6 months.
Note	For research use only
Clonality	Monoclonal
Source	Mouse Monoclonal Ab
Uniprot ID	Q14683
Dilution Range	Working with enzyme-conjugated streptavidin, 1: 50 - 1: 1,000
Expiration Date	12 months from date of receipt.