

Product Datasheet

HDAC4/HDAC5/HDAC9 (Ab-246/259/220) Conjugated Antibody (orb1629204)

Description	HDAC4/HDAC5/HDAC9 (Ab-246/259/220) Conjugated Antibody
Species/Host	Rabbit
Conjugation	Biotin
Storage	Store at 4°C for 6 months.
Note	For research use only
Clonality	Recombinant
Source	Rabbit Polyclonal Ab
Uniprot ID	Q9UQL6 , Q9UKV0 , P56524
Dilution Range	Working with enzyme-conjugated streptavidin, 1: 50 - 1: 1,000
Expiration Date	12 months from date of receipt.