

Product Datasheet

cofilin (Ab-3) Conjugated Antibody (orb1629346)

Description cofilin (Ab-3) Conjugated Antibody

Species/Host Rabbit

Conjugation Biotin

Storage Store at 4°C for 6 months.

Note For research use only

Clonality Recombinant

Source Rabbit Polyclonal Ab

Uniprot ID **P23528**

Dilution Range Working with enzyme-conjugated streptavidin, 1: 50 - 1: 1,000

Expiration Date 12 months from date of receipt.