

## Product Datasheet

### Akt2 (Ab-474) Conjugated Antibody (orb1629350)

**Description** Akt2 (Ab-474) Conjugated Antibody

**Species/Host** Rabbit

**Conjugation** Biotin

**Storage** Store at 4°C for 6 months.

**Note** For research use only

**Clonality** Recombinant

**Source** Rabbit Polyclonal Ab

**Uniprot ID** [P31751](#)

**Dilution Range** Working with enzyme-conjugated streptavidin, 1: 50 - 1: 1,000

**Expiration Date** 12 months from date of receipt.