

## Product Datasheet

### ASK1 (Ab-83) Conjugated Antibody (orb1629364)

<b>Description</b>	ASK1 (Ab-83) Conjugated Antibody
<b>Species/Host</b>	Rabbit
<b>Conjugation</b>	Biotin
<b>Storage</b>	Store at 4°C for 6 months.
<b>Note</b>	For research use only
<b>Clonality</b>	Recombinant
<b>Source</b>	Rabbit Polyclonal Ab
<b>Uniprot ID</b>	<a href="#">Q99683</a>
<b>Dilution Range</b>	Working with enzyme-conjugated streptavidin, 1: 50 - 1: 1,000
<b>Expiration Date</b>	12 months from date of receipt.