

Product Datasheet

PDK1 (Ab-241) Conjugated Antibody (orb1629403)

Description PDK1 (Ab-241) Conjugated Antibody

Species/Host Rabbit

Conjugation Biotin

Storage Store at 4°C for 6 months.

Note For research use only

Clonality Recombinant

Source Rabbit Polyclonal Ab

Uniprot ID [O15530](#)

Dilution Range Working with enzyme-conjugated streptavidin, 1: 50 - 1: 1,000

Expiration Date 12 months from date of receipt.