

Product Datasheet

Ah Receptor (phospho-Ser36) Conjugated Antibody (orb1629427)

| | |
|------------------------|---|
| Description | Ah Receptor (phospho-Ser36) Conjugated Antibody |
| Species/Host | Rabbit |
| Conjugation | Biotin |
| Storage | Store at 4°C for 6 months. |
| Note | For research use only |
| Clonality | Recombinant |
| Source | Rabbit Polyclonal Ab |
| Uniprot ID | P35869, A9YTQ3 |
| Dilution Range | Working with enzyme-conjugated streptavidin, 1: 50 - 1: 1,000 |
| Expiration Date | 12 months from date of receipt. |