

Product Datasheet

Histone H2A (Phospho-Thr121) Polyclonal Polyclonal Conjugated Antibody (orb1630337)

Description	Histone H2A (Phospho-Thr121) Polyclonal Polyclonal Conjugated Antibody
Species/Host	Rabbit
Conjugation	Biotin
Storage	Store at 4°C for 6 months.
Note	For research use only
Clonality	Recombinant
Source	Rabbit Polyclonal Ab
Uniprot ID	Q96KK5 , Q7L7L0 , P0C0S8 , P04908 , P20671 , Q6FI13
Dilution Range	Working with enzyme- conjugated streptavidin, 1: 50 - 1: 1,000
Expiration Date	12 months from date of receipt.