

Product Datasheet

Flt-1 (Phospho-Tyr1048) Polyclonal Polyclonal Conjugated Antibody (orb1630338)

| | |
|------------------------|---|
| Description | Flt-1 (Phospho-Tyr1048) Polyclonal Polyclonal Conjugated Antibody |
| Species/Host | Rabbit |
| Conjugation | Biotin |
| Storage | Store at 4°C for 6 months. |
| Note | For research use only |
| Clonality | Recombinant |
| Source | Rabbit Polyclonal Ab |
| Uniprot ID | P17948 |
| Dilution Range | Working with enzyme- conjugated streptavidin, 1: 50 - 1: 1,000 |
| Expiration Date | 12 months from date of receipt. |