

Product Datasheet

DNAM-1 (Phospho-Ser329) Polyclonal Polyclonal Conjugated Antibody (orb1630352)

Description	DNAM-1 (Phospho-Ser329) Polyclonal Polyclonal Conjugated Antibody
Species/Host	Rabbit
Conjugation	Biotin
Storage	Store at 4°C for 6 months.
Note	For research use only
Clonality	Recombinant
Source	Rabbit Polyclonal Ab
Uniprot ID	Q15762
Dilution Range	Working with enzyme- conjugated streptavidin, 1: 50 - 1: 1,000
Expiration Date	12 months from date of receipt.