

Product Datasheet

MerTK/Tyro3 (Phospho-Tyr749/681) Polyclonal Polyclonal Conjugated Antibody (orb1630369)

Description	MerTK/Tyro3 (Phospho-Tyr749/681) Polyclonal Polyclonal Conjugated Antibody
Species/Host	Rabbit
Conjugation	Biotin
Storage	Store at 4°C for 6 months.
Note	For research use only
Clonality	Recombinant
Source	Rabbit Polyclonal Ab
Uniprot ID	Q06418 , Q12866
Dilution Range	Working with enzyme-conjugated streptavidin, 1: 50 - 1: 1,000
Expiration Date	12 months from date of receipt.