

Product Datasheet

EphB1/2 (Phospho-Tyr594/604) Polyclonal Polyclonal Conjugated Antibody (orb1630412)

Description	EphB1/2 (Phospho-Tyr594/604) Polyclonal Polyclonal Conjugated Antibody
Species/Host	Rabbit
Conjugation	Biotin
Storage	Store at 4°C for 6 months.
Note	For research use only
Clonality	Recombinant
Source	Rabbit Polyclonal Ab
Uniprot ID	P54762, P29323
Dilution Range	Working with enzyme- conjugated streptavidin, 1: 50 - 1: 1,000
Expiration Date	12 months from date of receipt.