

Product Datasheet

Intestinal Cell Kinase (Phospho-Tyr159) Polyclonal Polyclonal Conjugated Antibody (orb1630428)

Description	Intestinal Cell Kinase (Phospho-Tyr159) Polyclonal Polyclonal Conjugated Antibody
Species/Host	Rabbit
Conjugation	Biotin
Storage	Store at 4°C for 6 months.
Note	For research use only
Clonality	Recombinant
Source	Rabbit Polyclonal Ab
Uniprot ID	Q9UPZ9
Dilution Range	Working with enzyme-conjugated streptavidin, 1: 50 - 1: 1,000
Expiration Date	12 months from date of receipt.