

## Product Datasheet

# RSK1/2/3/4 (Phospho-Ser221/227/S218/232) Polyclonal Conjugated Antibody (orb1630797)

|                        |   |
|------------------------|---|
| <b>Description</b>     | RSK1/2/3/4 (Phospho-Ser221/227/S218/232) Polyclonal Conjugated Antibody                           |
| <b>Species/Host</b>    | Rabbit  |
| <b>Conjugation</b>     | Biotin  |
| <b>Storage</b>         | Store at 4°C for 6 months.  |
| <b>Note</b>            | For research use only   |
| <b>Clonality</b>       | Recombinant   |
| <b>Source</b>          | Rabbit Polyclonal Ab  |
| <b>Uniprot ID</b>      | <a href="#">P51812</a> , <a href="#">Q9UK32</a> , <a href="#">Q15418</a> , <a href="#">Q15349</a> |
| <b>Dilution Range</b>  | Working with enzyme-conjugated streptavidin, 1: 50 - 1: 1,000                                     |
| <b>Expiration Date</b> | 12 months from date of receipt.   |