

## Product Datasheet

# GATA2 (Phospho-Ser401) Polyclonal Conjugated Antibody (orb1630844)

|                        |                                                                      |
|------------------------|----------------------------------------------------------------------|
| <b>Description</b>     | GATA2 (Phospho-Ser401)<br>Polyclonal Conjugated Antibody             |
| <b>Species/Host</b>    | Rabbit                                                               |
| <b>Conjugation</b>     | Biotin                                                               |
| <b>Storage</b>         | Store at 4°C for 6 months.                                           |
| <b>Note</b>            | For research use only                                                |
| <b>Clonality</b>       | Recombinant                                                          |
| <b>Source</b>          | Rabbit Polyclonal Ab                                                 |
| <b>Uniprot ID</b>      | <b>P23769</b>                                                        |
| <b>Dilution Range</b>  | Working with enzyme-<br>conjugated streptavidin, 1: 50 -<br>1: 1,000 |
| <b>Expiration Date</b> | 12 months from date of<br>receipt.                                   |