

Product Datasheet

FOXO1/3/4-pan (Phospho-Thr24/32) Polyclonal Conjugated Antibody (orb1630863)

Description	FOXO1/3/4-pan (Phospho-Thr24/32) Polyclonal Conjugated Antibody
Species/Host	Rabbit
Conjugation	Biotin
Storage	Store at 4°C for 6 months.
Note	For research use only
Clonality	Recombinant
Source	Rabbit Polyclonal Ab
Uniprot ID	Q12778 , O43524 , P98177
Dilution Range	Working with enzyme-conjugated streptavidin, 1: 50 - 1: 1,000
Expiration Date	12 months from date of receipt.