

Product Datasheet

Acetyl-CoA Carboxylase (Phospho-Ser79) Polyclonal Conjugated Antibody (orb1630887)

Description	Acetyl-CoA Carboxylase (Phospho-Ser79) Polyclonal Conjugated Antibody
Species/Host	Rabbit
Conjugation	Biotin
Storage	Store at 4°C for 6 months.
Note	For research use only
Clonality	Recombinant
Source	Rabbit Polyclonal Ab
Uniprot ID	Q13085
Dilution Range	Working with enzyme- conjugated streptavidin, 1: 50 - 1: 1,000
Expiration Date	12 months from date of receipt.