

Product Datasheet

EGFR (phospho-Tyr1197) Goat Polyclonal Polyclonal Conjugated Antibody (orb1630903)

| | |
|------------------------|---|
| Description | EGFR (phospho-Tyr1197) Goat Polyclonal Polyclonal Conjugated Antibody |
| Species/Host | Goat |
| Conjugation | Biotin |
| Storage | Store at 4°C for 6 months. |
| Note | For research use only |
| Clonality | Polyclonal |
| Source | Goat Polyclonal Ab |
| Uniprot ID | P00533 |
| Dilution Range | Working with enzyme- conjugated streptavidin, 1: 50 - 1: 1,000 |
| Expiration Date | 12 months from date of receipt. |