

## Product Datasheet

# Integrin $\beta$ 3 (Phospho-Tyr785) Polyclonal Conjugated Antibody (orb1630959)

<b>Description</b>	Integrin $\beta$ 3 (Phospho-Tyr785) Polyclonal Conjugated Antibody
<b>Species/Host</b>	Rabbit
<b>Conjugation</b>	Biotin
<b>Storage</b>	Store at 4°C for 6 months.
<b>Note</b>	For research use only
<b>Clonality</b>	Recombinant
<b>Source</b>	Rabbit Polyclonal Ab
<b>Uniprot ID</b>	<b>P05106</b>
<b>Dilution Range</b>	Working with enzyme- conjugated streptavidin, 1: 50 - 1: 1,000
<b>Expiration Date</b>	12 months from date of receipt.