

## Product Datasheet

# CREB (Phospho-Ser129) Polyclonal Conjugated Antibody (orb1630964)

|                        |  |
|------------------------|--|
| <b>Description</b>     | CREB (Phospho-Ser129)<br>Polyclonal Conjugated Antibody              |
| <b>Species/Host</b>    | Rabbit   |
| <b>Conjugation</b>     | Biotin   |
| <b>Storage</b>         | Store at 4°C for 6 months.   |
| <b>Note</b>            | For research use only  |
| <b>Clonality</b>       | Recombinant  |
| <b>Source</b>          | Rabbit Polyclonal Ab   |
| <b>Uniprot ID</b>      | <b>P16220</b>  |
| <b>Dilution Range</b>  | Working with enzyme-<br>conjugated streptavidin, 1: 50 -<br>1: 1,000 |
| <b>Expiration Date</b> | 12 months from date of receipt.                                      |