

Product Datasheet

SAPK/JNK (Phospho-Thr183) Polyclonal Conjugated Antibody (orb1630980)

Description	SAPK/JNK (Phospho-Thr183) Polyclonal Conjugated Antibody
Species/Host	Rabbit
Conjugation	Biotin
Storage	Store at 4°C for 6 months.
Note	For research use only
Clonality	Recombinant
Source	Rabbit Polyclonal Ab
Uniprot ID	P53779 , P45984
Dilution Range	Working with enzyme- conjugated streptavidin, 1: 50 - 1: 1,000
Expiration Date	12 months from date of receipt.