

Product Datasheet

Synaptotagmin 1/2 (Phospho-Thr202/199) Polyclonal Conjugated Antibody (orb1631005)

Description	Synaptotagmin 1/2 (Phospho-Thr202/199) Polyclonal Conjugated Antibody
Species/Host	Rabbit
Conjugation	Biotin
Storage	Store at 4°C for 6 months.
Note	For research use only
Clonality	Recombinant
Source	Rabbit Polyclonal Ab
Uniprot ID	Q8N9I0 , P21579
Dilution Range	Working with enzyme-conjugated streptavidin, 1: 50 - 1: 1,000
Expiration Date	12 months from date of receipt.