

## Product Datasheet

### General Testosterone ELISA Kit (orb1631167)

<b>Description</b>	General Testosterone ELISA Kit
<b>Reactivity</b>	All
<b>Tested Applications</b>	ELISA
<b>Note</b>	For research use only
<b>Application notes</b>	Detect Range: 3.9--250 pg/mL. Assay Time: 3-5 h. Detection wavelength: 450 nm. Sample Volume: 50 T. Sample Type: Serum, plasma, tissue homogenates and other biological fluids.
<b>Expiration Date</b>	Please enquire.