

## Product Datasheet

### Recombinant Rat FGF-basic (orb1631602)

<b>Description</b>	Recombinant Rat FGF-basic
<b>Form/Appearance</b>	Lyophilized from a 0.2 µm filtered solution in PBS, pH 7.4.
<b>Storage</b>	This lyophilized preparation is stable at 2-8 °C, but should be kept at -20 °C for long term storage, preferably desiccated. Upon reconstitution, the preparation is stable for up to one week at 2-8 °C. For maximal stability, apportion the reconstituted preparation into working aliquots and store at -20 °C to -70 °C. Avoid repeated freeze/thaw cycles.
<b>Note</b>	For research use only
<b>Purity</b>	> 98 % by SDS-PAGE and HPLC analyses.
<b>Source</b>	Rat
<b>MW</b>	Approximately 16.4 kDa, a single non-glycosylated polypeptide chain containing 146 amino acids.
<b>Uniprot ID</b>	<b>P13109</b>
<b>Expiration Date</b>	6 months from date of receipt.