

Product Datasheet

Mouse Pcsk9 protein (orb163426)

Description	Recombinant Mouse Proprotein Convertase Subtilisin/Kexin Type 9 is produced by our Mammalian express...
Endotoxins	Less than 0.1 ng/ug (1 EU/ug)
Form/Appearance	Supplied as a 0.2 um filtered solution of 50mM HEPES, 150mM NaCl, 20% Glycerol, pH 7.4.
Storage	Store at -20°C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.
Tag	C-His
Note	For research use only
Application notes	6His tag at the C-terminus.
Uniprot ID	Q80W65
Expiration Date	6 months from date of receipt.