

## Product Datasheet

# SARS-CoV-2 (2019-nCoV) Spike Protein (S2 ECD, His tag) (orb1635314)

<b>Description</b>	SARS-CoV-2 (2019-nCoV) Spike Protein (S2 ECD, His tag)
<b>Species/Host</b>	Virus
<b>Storage</b>	Samples are stable for up to twelve months from date of receipt at -20°C to -80°C Store it under sterile conditions at -20°C to -80°C. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.
<b>Note</b>	For research use only
<b>Expiration Date</b>	6 months from date of receipt.