

Product Datasheet

Human RHEB protein (orb163538)

Description	Recombinant Human Ras homolog enriched in brain is produced by our E.coli expression system and the ...
Endotoxins	Less than 0.1 ng/ug (1 EU/ug)
Form/Appearance	Lyophilized from a 0.2 um filtered solution of 20mM Tris, 10mM GSH, pH 8.0.
Storage	Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7°C for 2-7 days. Aliquots of reconstituted samples are stable at -20°C for 3 months.
Tag	N-GST
Note	For research use only
Application notes	GST tag at the N-terminus.
Uniprot ID	Q15382
Expiration Date	6 months from date of receipt.