

## Product Datasheet

# Dolichyl-diphosphooligosaccharide--protein glycosyltransferase 48 kDa subunit (orb1637903)

<b>Description</b>	Dolichyl-diphosphooligosaccharide--protein glycosyltransferase 48 kDa subunit
<b>Species/Host</b>	E. coli
<b>Storage</b>	Store at -20°C. (Avoid repeated freezing and thawing.)
<b>Tag</b>	His
<b>Note</b>	For research use only
<b>Expiration Date</b>	6 months from date of receipt.