

## Product Datasheet

# Phosphatidylinositol-binding clathrin assembly protein (orb1637955)

|                        |  |
|------------------------|--|
| <b>Description</b>     | Phosphatidylinositol-binding clathrin assembly protein |
| <b>Species/Host</b>    | E. coli  |
| <b>Storage</b>         | Store at -20°C. (Avoid repeated freezing and thawing.) |
| <b>Tag</b>             | His  |
| <b>Note</b>            | For research use only                                  |
| <b>Expiration Date</b> | 6 months from date of receipt.                         |