

Product Datasheet

SEMA4A antibody (orb164341)

Description	Rabbit polyclonal antibody to SEMA4A (N-term).
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	SEMA4A
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	83574
Uniprot ID	Q9H3S1
NCBI	NP_001180229.1 , NP_001180231.1 , NP_071762.2 , NP_001180230.1
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.