

Product Datasheet

Recombinant Escherichia coli Iron-sulfur cluster repair protein ytfE (ytfE) (Active) (orb1646259)

Description	Recombinant Escherichia coli Iron-sulfur cluster repair protein ytfE (ytfE) (Active)
Species/Host	E. coli
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Tag	N-terminal GST-tagged
Note	For research use only
Expiration Date	6 months from date of receipt.