

## Product Datasheet

# Recombinant Human Cysteine-rich secretory protein LCCL domain-containing 2(CRISPLD2) (orb1646443)

<b>Description</b>	Recombinant Human Cysteine-rich secretory protein LCCL domain-containing 2(CRISPLD2)
<b>Species/Host</b>	Virus
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Tag</b>	N-terminal 10xHis-tagged and C-terminal Myc-tagged
<b>Note</b>	For research use only
<b>Expiration Date</b>	6 months from date of receipt.