

Product Datasheet

Recombinant Human Cysteine-rich secretory protein LCCL domain-containing 2(CRISPLD2) (orb1646451)

Description	Recombinant Human Cysteine-rich secretory protein LCCL domain-containing 2(CRISPLD2)
Species/Host	Virus
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Tag	N-terminal MBP-tagged and C-terminal 6xHis-tagged
Note	For research use only
Expiration Date	6 months from date of receipt.