

Product Datasheet

Recombinant Mycobacterium tuberculosis Immunogenic protein MPT64(mpt64) (orb1646479)

| | |
|------------------------|---|
| Description | Recombinant Mycobacterium tuberculosis Immunogenic protein MPT64(mpt64) |
| Species/Host | Virus |
| Storage | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |
| Tag | N-terminal 10xHis-tagged |
| Note | For research use only |
| Expiration Date | 6 months from date of receipt. |