

Product Datasheet

Recombinant Mouse Disintegrin and metalloproteinase domain-containing protein 12(Adam12), partial (orb1646500)

Description	Recombinant Mouse Disintegrin and metalloproteinase domain-containing protein 12(Adam12), partial
Species/Host	Virus
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Tag	N-terminal 10xHis-tagged and C-terminal Myc-tagged
Note	For research use only
Expiration Date	6 months from date of receipt.