

## Product Datasheet

# Recombinant Human Receptor-interacting serine/threonine-protein kinase 1(RIPK1),partial (orb1646551)

<b>Description</b>	Recombinant Human Receptor-interacting serine/threonine-protein kinase 1(RIPK1),partial
<b>Species/Host</b>	E. coli
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Tag</b>	N-terminal 10xHis-tagged and C-terminal Myc-tagged
<b>Note</b>	For research use only
<b>Expiration Date</b>	6 months from date of receipt.