

Product Datasheet

Recombinant Mouse Endothelin-converting enzyme-like 1(Ecel1),partial (orb1646662)

Description	Recombinant Mouse Endothelin-converting enzyme-like 1(Ecel1),partial
Species/Host	E. coli
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Tag	N-terminal 10xHis-tagged and C-terminal Myc-tagged
Note	For research use only
Expiration Date	6 months from date of receipt.