

Product Datasheet

TPM2 antibody (orb164752)

Description	Rabbit polyclonal antibody to TPM2 (C-term).
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	TPM2
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	32851
Uniprot ID	P07951
NCBI	NP_003280.2 , NP_998839.1
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.