

## Product Datasheet

# Recombinant Human immunodeficiency virus type 2 subtype A Protein Vpx(vpx) (orb1647826)

<b>Description</b>	Recombinant Human immunodeficiency virus type 2 subtype A Protein Vpx(vpx)
<b>Species/Host</b>	Yeast
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Tag</b>	N-terminal 6xHis-tagged
<b>Note</b>	For research use only
<b>Expiration Date</b>	6 months from date of receipt.