

## Product Datasheet

# Recombinant Human V-Set and Ig Domain-Containing Protein 8/VSIG8 (C-6His) (orb1648085)

<b>Description</b>	Recombinant Human V-Set and Ig Domain-Containing Protein 8/VSIG8 (C-6His)
<b>Species/Host</b>	Human
<b>Storage</b>	Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7°C for 2-7 days. Aliquots of reconstituted samples are stable at -20°C for 3 months.
<b>Note</b>	For research use only
<b>Expiration Date</b>	6 months from date of receipt.