

Product Datasheet

Recombinant Mouse Interleukin-18 Binding Protein Isoform d/IL-18 BPd (C-Fc) (orb1648185)

Description	Recombinant Mouse Interleukin-18 Binding Protein Isoform d/IL-18 BPd (C-Fc)
Species/Host	Human
Storage	Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7°C for 2-7 days. Aliquots of reconstituted samples are stable at -20°C for 3 months.
Note	For research use only
Expiration Date	6 months from date of receipt.