

Product Datasheet

Recombinant Human/Mouse/Rat Irisin/FNDC5 (C-6His) (orb1648194)

Description	Recombinant Human/Mouse/Rat Irisin/FNDC5 (C-6His)
Species/Host	Human
Storage	Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7°C for 2-7 days. Aliquots of reconstituted samples are stable at -20°C for 3 months.
Note	For research use only
Expiration Date	6 months from date of receipt.