

## Product Datasheet

# Recombinant Human Pro-Neuregulin-1/NRG1- $\beta$ 1/HRG1- $\beta$ 1 (Met1-Lys246, N-6His) (orb1648455)

<b>Description</b>	Recombinant Human Pro-Neuregulin-1/NRG1- $\beta$ 1/HRG1- $\beta$ 1 (Met1-Lys246, N-6His)
<b>Species/Host</b>	E. coli
<b>Storage</b>	Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7°C for 2-7 days. Aliquots of reconstituted samples are stable at -20°C for 3 months.
<b>Note</b>	For research use only
<b>Expiration Date</b>	6 months from date of receipt.