

Product Datasheet

Recombinant Human Pro-Neuregulin-1/NRG1- β 1/HRG1- β 1 (Thr176-Lys246) (orb1648456)

Description	Recombinant Human Pro-Neuregulin-1/NRG1- β 1/HRG1- β 1 (Thr176-Lys246)
Species/Host	E. coli
Storage	Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7°C for 2-7 days. Aliquots of reconstituted samples are stable at -20°C for 3 months.
Note	For research use only
Expiration Date	6 months from date of receipt.